

Patient Information

Patient ID: [REDACTED]

Birthdate: [REDACTED]

Test started: August 11, 2023, 1:53 PM

Report generated: August 11, 2023, 2:06 PM

Device and Test Information

RETeval™


Serial number: [REDACTED]

Test protocol: DR Assessment


Manufacturer: LKC Technologies, Inc.

Firmware version: 2.13.5 Reference data: 2023.23 6966f91


Electrodes: Sensor Strips



DR Score	19.8
Operator-selected limits (7.0 ↔ 23.4)	Within limits
95% Reference interval (14.3 ↔ 22.6)	68%




Right Eye




Left Eye

ERG

	ms		μV			ms		μV	
16 Td·s	26.8 (16%)	25.9 ↔ 31.7	21.1 (40%)	12.3 ↔ 38.4	16 Td·s	26.6 (12%)	25.9 ↔ 31.7	21.3 (41%)	12.3 ↔ 38.4
32 Td·s	26.0 (11%)	25.2 ↔ 30.0	26.1 (45%)	15.3 ↔ 43.9	32 Td·s	25.9 (10%)	25.2 ↔ 30.0	27.2 (50%)	15.3 ↔ 43.9




16 Td-s White 28.3 Hz flicker, no background



32 Td-s White 28.3 Hz flicker, no background


Pupil

Area ratio: 1.6 (16%) 1.4 ↔ 3.2



4.0 Td-s White 28.3 Hz flicker, no background

Area ratio: 1.5 (7%) 1.4 ↔ 3.2



32 Td-s White 28.3 Hz flicker, no background